


BREAKTHROUGH ANALYSIS ON THE AUTOCHEM



INTRODUCTION

The AutoChem is a flow-through system that allows for the characterization of catalysts for surface activity. The AutoChem can also function as a basic single component breakthrough unit allowing for carrier, preparation, and analysis gas flows. This document will describe how to perform a breakthrough analysis using the AutoChem alongside a mass spectrometer.

RESULTS

The results of the breakthrough experiment are shown below. Zeolite 13X showed strong adsorption for CO₂ at ambient conditions reaching a capacity of 4.53 mmol/g. Additionally, the breakthrough curve is steep signifying that there are little to no mass transfer limitations in the system.




Figure 2. CO₂ breakthrough adsorption in zeolite 13X.

CONCLUSIONS

The AutoChem is a capable unit for performing simple single component breakthrough analysis. This was confirmed via analysis of zeolite 13X, achieving a CO₂ adsorption capacity of 4.53 mmol/g at 30 °C and a flowrate of 20 sccm.