

# ***GEOPYC® 1365***

---

*ENVELOPE DENSITY ANALYZER*



***CALCULATIONS***

---

Nov 2016  
(Ver-)

## *Table of Contents*

Calculated Blank Counts _____	1
Calculated Volume Calibration _____	1
Force Calibration _____	2
Percent of Sample Volume _____	2
Percent Porosity _____	3
Plunger Position _____	3
Sample Envelope Density _____	4
Specific Pore Volume _____	4
Volume Calibration _____	5
Zero Depth Position _____	5
Compendium of Variables _____	6

## ***CALCULATED BLANK COUNTS***

---

Assume that multiple blank data sets are available for a given chamber. Each set is required to have a unique chamber+medium mass. When an analysis using calculated blank counts is started, the blank data set index  $j$  is found such that  $m_j \leq m < m_{j+1}$  where  $m$  is the chamber+medium mass for the sample to be analyzed. An error is signaled if there is no such index.

For each measurement cycle to be done the calculated blank counts is found by linearly interpolating counts in the corresponding cycles of blank sets  $j$  and  $j+1$  as a function of chamber+medium mass.

$$n_i = n_{i,j} + \left( n_{i,j+1} - n_{i,j} \right) \frac{m - m_j}{m_{j+1} - m_j}$$

Each set may have a different number of cycles,  $N$ , but an error is signaled if

$$N > \min(N_j, N_{j+1})$$

## ***CALCULATED VOLUME CALIBRATION***

---

A scale factor can be calculated from the geometry of the chamber.

$$\begin{aligned} k_V &= \Delta V / \Delta x \\ &= A \\ &= \pi d^2 / 4 \end{aligned}$$

We need to divide by 10 to express the scale factor in  $\text{cm}^3/\text{mm}$ .

$$k_V = \pi d^2 / 40$$

## FORCE CALIBRATION

---

The force transducer is calibrated using a spring with known spring constant.

First, the force signal and position counts are recorded at the home position ( $V_0$  at  $n_0$ ).

Plunger positions at two different measured forces are recorded,  $n_{low}$  and  $n_{high}$ . The plunger is then moved in 10 equal steps from  $n_{low}$  and  $n_{high}$ . Position,  $n_i$  and force signals  $V_i$  are measured at each step.

A least-squares straight-line fit is done through

$$x_i = V_i - V_0$$

$$y_i = k_s k_x (n_i - n_0)$$

where  $x_i$  is in volts, and  $y_i$  is in newtons. The new force calibration slope is the slope of the line in N/V.

## PERCENT OF SAMPLE VOLUME

---

The fraction of the volume of the contents of the chamber that is due to the sample is

$$f_s = V_s / (V_s + V_b)$$

where

$$\begin{aligned} V_b &= \text{volume of the chamber contents without sample} \\ V_s &= \text{volume of the chamber contents with sample} \end{aligned}$$

The calibration scale factor divides out so we can write

$$f_s = (n_b - n_s) / (n_0 - n_s)$$

The zero depth position,  $n_0$ , is needed since the total volume is relative to the empty chamber. The percent of sample volume reported is  $100\% \bar{f}_s$  where  $\bar{f}_s$  is calculated from the average of  $n_b$  and  $n_s$  over the consolidation cycles.

## ***PERCENT POROSITY***

---

Porosity is determined from envelope density and entered absolute density.

$$P = 100\% \left(1 - \rho_e / \rho_a\right)$$

## ***PLUNGER POSITION***

---

The plunger is moved by a screw attached to a stepper motor. Distance moved in millimeters is proportional to the number of counts from the stepper motor encoder.

$$x = k_x n$$

$$r_m = 400 \text{ counts/revolution}$$

The lead screw is 5/16-24 UNF.

$$r_t = 24 \text{ revolution/in}$$

$$k_x = (r_m r_t)^{-1} \text{ in/count}$$

$$k_x = 1 / 9600 = \text{in/count} = 127 / 4800 = \text{mm/count}$$

## **SAMPLE ENVELOPE DENSITY**

---

Envelope volume is the volume of an object that includes small pores and cavities. The GeoPyc reports volume that is not accessible to the Dry Flo medium as envelope volume.

When compressed by the piston, the volume of material in the chamber is  $V_e = Ax$  where  $A$  is the cross-sectional area of the chamber and  $x$  is the plunger position relative to the bottom of the chamber.

Sample volume is found as the difference in volume of the sample and medium and the volume of a blank run.

$$\begin{aligned} V_e &= A((x_b - x_0) - (x_s - x_0)) \\ &= A(x_b - x_s) \\ &= Ak_x(n_b - n_s) \end{aligned}$$

This voids the need to know the position of the bottom of the chamber.

In practice, a scale factor,  $k_V$ , determined by volume calibration is used instead of  $Ak_x$ .

$$V_e = k_V(n_b - n_s)$$

Envelope density is found from the measured envelope volume and the entered sample mass.

$$\rho_e = m / V_e$$

## **SPECIFIC PORE VOLUME**

---

Specific pore volume is determined from envelope density and entered absolute density.

$$V_p = (1 / \rho_e - 1 / \rho_a)$$

## VOLUME CALIBRATION

---

Measurement of a reference sample with a known properties allows  $k_V$  to be determined.

If the reference volume,  $V_r$ , is known the envelope volume equation can be solved for  $k_V$ . The average over the measurement cycles is used.

$$k_V = V_r / (\bar{n}_b - \bar{n}_s)$$

This can be written in terms of the average envelope volume from the calibration run.

$$k_V = k_V^* \frac{V_r}{V_e}$$

where  $k_V^*$  is the conversion factor used in the run.

If mass, absolute density, and porosity are known, the reference volume can be calculated from the porosity equation.

$$\begin{aligned} P &= 1 - \rho_e / \rho_a \\ &= 1 - \frac{m}{V_e \rho_a} \end{aligned}$$

$$V_r = \frac{m}{\rho_a(1 - P)}$$

Volume error is the difference between the average measured envelope volume and the entered or calculated reference volume.

$$\Delta V = V_e - V_r$$

## ZERO DEPTH POSITION

---

When finding the plunger position at the bottom of the chamber, the plunger is assumed to be extended by 1/2" in (4800 counts).

$$n_0 = n + 4800$$

---

## COMPENDIUM OF VARIABLES

---

$k_s$	=	spring constant
$k_V$	=	scale factor for converting plunger position to volume
$k_x$	=	scale factor for converting plunger position in counts to millimeters
$m$	=	sample mass
$n$	=	plunger position in stepper motor counts
$n_0$	=	plunger counts at the bottom of the chamber
$n_b$	=	plunger counts without sample
$n_s$	=	plunger counts with sample
$P$	=	percent porosity
$V_e$	=	envelope volume
$V_p$	=	pore volume
$V_r$	=	calibration reference volume
$x$	=	plunger position in millimeters
$x_0$	=	plunger position at the bottom of the chamber
$x_b$	=	plunger position without sample
$x_s$	=	plunger position with sample
$\rho_a$	=	absolute density
$\rho_e$	=	envelope density